

FALCON™

Anterior Cervical Plate System



FALCON™



Anterior Cervical Plate System

The Falcon™ Cervical Plate System designed to simplify anterior cervical spinal fusion with a no step locking mechanism. The Falcon cervical plate features a sleek design for easy insertion and large graft windows for enhanced visual confirmation.


Features:

- Large window for optimal visualization
- Simple no -step locking clip for fast screw insertion and locking
- Self-drilling and self-tapping, fixed and variable angle screws for maximum purchase
- Large assortment of plate lengths to accommodate a variety of patient anatomies (1-4 level Plates)
- Screw lengths: 10-16mm Variable and Fixed Angle screws, self-drilling and self-tapping

Variable Screw



Falcon Plate Sizes



Fixed Screw



Self Drilling & Self Tapping Screws



400 Erin Drive, Knoxville, TN 37919 | O: 865.246.3333 | F: 865.246.3334 | choicespine.com