

FALCON™

CERVICAL PLATE SYSTEM

FAST, REPRODUCIBLE RESULTS



FALCON™

The FALCON™ Cervical Plating System simplifies anterior cervical spinal fusion with an improved semi-constrained system that offers fast, reproducible results. The Falcon cervical plate features a sleek, contemporary design for easy insertion and large graft windows for enhanced visual confirmation of the allograft. The instrumentation facilitates complete control by the surgeon during the procedure.

KEY PLATE FEATURES

- Large Graft Window
- Simple, One-Step Locking Mechanism
- Low Profile Design
- Self-Drilling, Self-Tapping, Fixed & Variable Angle Screws
- Large Assortment of Plate Lengths
- Scalloped Sides for Lateral Visualization



4 Level 80mm - 84mm (4mm inc.)



3 Level 40mm - 62mm (4mm inc.)



2 Level 20mm - 44mm (2mm inc.)



1 Level 10mm - 26mm (2mm inc.)



Item #	Implant Description	Available Lengths
FT10-BR410, FT10-BR412, FT10-BR414, FT10-BR416	4.0mm Fixed Angle Self Drilling Screws	10mm, 12mm, 14mm, 16mm
FT10-BP410, FT10-BP412, FT10-BP414, FT10-BP416	4.0mm Variable Angle Self Drilling Screws	10mm, 12mm, 14mm, 16mm
FT10-BN410, FT10-BN412, FT10-BN414, FT10-BN416	4.5mm Fixed Angle Self Tapping Screws	10mm, 12mm, 14mm, 16mm
FT10-BL410, FT10-BL412, FT10-BL414, FT10-BL416	4.5mm Variable Angle Self Tapping Screws	10mm, 12mm, 14mm, 16mm
FT10-BH410, FT10-BH412, FT10-BH414, FT10-BH416	4.0mm Fixed Angle Self Tapping Screws	10mm, 12mm, 14mm, 16mm
FT10-BF410, FT10-BF412, FT10-BF414, FT10-BF416	4.0mm Variable Angle Self Tapping Screws	10mm, 12mm, 14mm, 16mm
FT10-BE412, FT10-BE414, FT10-BE416	3.7mm Variable Angle Self Drilling Screws	12mm, 14mm, 16mm

ChoiceSpine™
Propelling Spinal Surgery

400 Erin Drive
Knoxville, TN 37919
865.246.3333 office 865.246.3334 fax
www.ChoiceSpine.com

Lit # Falcon SS
REV 02
7/17