

SHARK™ Lumbar Spacer



Features / Benefits

- Self-distracting nose for ease of insertion
- Asymmetric teeth for optimizing endplate engagement
- Large graft window for more bone graft volume

SHARK™

Lumbar Spacer

Parallel (0°)

22mm Length		25mm Length		28mm Length		32mm Length	
Heights	Catalog No.	Heights	Catalog No.	Heights	Catalog No.	Heights	Catalog No.
8mm	PP20-22080	8mm	PP20-25080	8mm	PP20-28080	8mm	PP20-32080
9mm	PP20-22090	9mm	PP20-25090	9mm	PP20-28090	9mm	PP20-32090
10mm	PP20-22100	10mm	PP20-25100	10mm	PP20-28100	10mm	PP20-32100
11mm	PP20-22110	11mm	PP20-25110	11mm	PP20-28110	11mm	PP20-32110
12mm	PP20-22120	12mm	PP20-25120	12mm	PP20-28120	12mm	PP20-32120
13mm	PP20-22130	13mm	PP20-25130	13mm	PP20-28130	13mm	PP20-32130
14mm	PP20-22140	14mm	PP20-25140	14mm	PP20-28140	14mm	PP20-32140

Lordotic (4°)

28mm Length	
Heights	Catalog No.
8mm	PP20-28084
9mm	PP20-28094
10mm	PP20-28104
11mm	PP20-28114
12mm	PP20-28124
13mm	PP20-28134
14mm	PP20-28144



400 Erin Drive
Knoxville, Tennessee 37919
865.246.3333 office | 865.246.3334 fax

Optional Sizes
Heights 7mm, 15mm-18mm Available Upon Request

www.ChoiceSpine.net

Please contact Choice Spine for a copy of this device's product information/
package insert for additional information, including indications for use.

Lit # SharkSS
Rev 01
09/13