




STEALTH™

PEEK Anterior
Cervical Spacer


STEALTH™

PEEK Anterior Cervical Spacer

The ChoiceSpine STEALTH PEEK Anterior Cervical Spacer implant is intended for anterior cervical spine intervertebral body fusion at one level from the C2-C3 disc space to the C7-T1 disc. The STEALTH spacer can be inserted in combination with the BOOMERANG™ Anterior Cervical Plate System for greater efficiency for our surgeon partners. Made from radiolucent PEEK, STEALTH comes in a variety of lengths, widths and heights to accommodate diverse patient anatomies.

FEATURES


- Beveled nose for easy insertion
- Prominent teeth for optimizing endplate engagement
- Large graft window for more bone graft volume
- 0° and 6° lordosis options
- 5, 6, 7, 8, 9, 10, 11, and 12 millimeter height options
- 12x14mm and 14x16mm footprints
- Compatible with Boomerang™ Anterior Cervical Plate



..... Boomerang™ Anterior Cervical Plate Compatible



..... Secure Threaded Inserter



..... Excellent Visibility Under Fluoroscopy

ChoiceSpine™
Propelling Spinal Surgery

400 Erin Drive | Knoxville, TN 37919
865.246.3333 | fax 865.246.3334
choicespine.com