


STRATOFUSE[®] CERVICAL

STRUCTURAL ALLOGRAFT



STRATOFUSE® CERVICAL

The STRATOFUSE Cervical Graft is 100% cortical bone. With anti-migration ribs and a wide footprint, this structural allograft maximizes disc space coverage while providing space to pack additional bone graft material. STRATOFUSE Cervical is saline packed to reduce graft fracture upon implantation and is available in various sizes to ensure clinical applicability.



FEATURES AND BENEFITS

- 100% cortical bone
- Large footprint: 13.4mm X 15.4mm
- Heights: 5mm – 10mm in 1mm increments
- 5° lordosis and parallel options
- Saline packed
- Anti-migration ribs
- Larger graft window than smaller footprint grafts
- Five year shelf life
- Sterile: SAL 10⁻⁶

Item #	Description	Size
BA65-131505	Cervical Graft, 5 Deg	5 mm
BA65-131506	Cervical Graft, 5 Deg	6 mm
BA65-131507	Cervical Graft, 5 Deg	7 mm
BA65-131508	Cervical Graft, 5 Deg	8 mm
BA65-131509	Cervical Graft, 5 Deg	9 mm
BA65-131510	Cervical Graft, 5 Deg	10mm
BA60-131505	Cervical Graft, Parallel	5 mm
BA60-131506	Cervical Graft, Parallel	6 mm
BA60-131507	Cervical Graft, Parallel	7 mm
BA60-131508	Cervical Graft, Parallel	8 mm
BA60-131509	Cervical Graft, Parallel	9 mm
BA60-131510	Cervical Graft, Parallel	10mm



A portion of all Biologics proceeds is donated to Smoky Mountain Service Dogs. This non-profit organization is dedicated to helping wounded veterans. Their mission is to enhance the physical and psychological quality of life for veterans with disabilities by providing custom-trained mobility-assistance service dogs. We are proud to partner with such an outstanding organization that gives back to those who served.



400 Erin Drive
Knoxville, TN 37919
865.246.3333 office 865.246.3334 fax
www.ChoiceSpine.com

Lit # Stratofuse Cervical SS
REV 02
09/18