


STRATOFUSE® DBM 100

100% DEMINERALIZED BONE MATRIX



STRATOFUSE® DBM 100

STRATOFUSE DBM 100 is a formable and flowable 100% demineralized bone matrix putty. STRATOFUSE DBM 100 does not contain any extrinsic carriers and is entirely derived from allograft bone. It resists irrigation and can be used to augment autologous cancellous or corticocancellous grafts and is provided in a ready-to-use syringe.

STRATOFUSE DBM 100 is a ready-to-use putty formulation that delivers excellent handling and timesaving convenience. Through proprietary methods, this allograft is processed to preserve naturally occurring BMPs and growth factors. This technique combined with a no-extrinsic-carrier formulation results in an end product with a higher concentration of endogenous growth factors than many DBMs available.



FEATURES AND BENEFITS

- Osteoinductive: each lot is tested in-vivo
- 100% Demineralized allograft bone, no extrinsic carrier
- Formable and irrigation resistant
- Osteoconductive
- Two year shelf life
- Sterile: SAL 10⁻⁶

Item #	Description	Size
BA20-000001	STRATOFUSE DBM 100,	1.0cc
BA20-000002	STRATOFUSE DBM 100,	2.5cc
BA20-000005	STRATOFUSE DBM 100,	5.0cc
BA20-000010	STRATOFUSE DBM 100,	10 cc



A portion of all Biologics proceeds is donated to Smoky Mountain Service Dogs. This non-profit organization is dedicated to helping wounded veterans. Their mission is to enhance the physical and psychological quality of life for veterans with disabilities by providing custom-trained mobility-assistance service dogs. We are proud to partner with such an outstanding organization that gives back to those who served.


ChoiceSpineTM
Propelling Spinal Surgery

400 Erin Drive
Knoxville, TN 37919
865.246.3333 office 865.246.3334 fax
www.ChoiceSpine.com

Lit # Stratofuse DBM 100 SS
REV 02
09/18