



STRATOFUSE® LUMBAR

Structural Lumbar Allograft





STRATOFUSE® LUMBAR

Structural Lumbar Allograft



STRATOFUSE® Lumbar Spacers are made of 100% cortical bone to maximize the strength and stability of each implant. Anti-migration ribs and a length of 23mm allows for placement of the allograft to provide the optimal disc space coverage. The multiple height options available and dual-surgery technique capability allows for physicians to match their clinical needs and desires.

Featured Benefits

- 100% Cortical Bone
- Length: 23mm, Depth: 9mm
- Heights: 7-15mm, 2mm Increments
- Anti-Migration Ribs
- Stored at Ambient Temperature
- Sterile: SAL 10^{-6}
- Source: Femur or Tibia

Available Footprints

- 7x23mm
- 9x23mm
- 11x23mm
- 13x23mm
- 15x23mm



Visit Biologics Portfolio

Spine the Right Way.SM

Distributed by ChoiceSpine LLC
400 Erin Drive
Knoxville, TN 37919
choicespine.com
865-243-3383
salessupport@choicespine.com

Copyright © 2022 ChoiceSpine LLC. All rights reserved.
Products referenced with ® are registered trademarks of ChoiceSpine LLC.

SS-000008 REV A 09/22