

TIGERSHARK™ THREADED TL

3D Printed Titanium Curved Interbody



TIGERSHARK™ THREADED TL

3D Printed Titanium Curved Interbody

ChoiceSpine's TigerShark™ Threaded TL is a 3D Printed curved interbody for TLIF procedures it is offered in many different footprints and 6° lordosis to accommodate patient anatomies. TigerShark Threaded TL features our propriety BioBond™, a porous structure technology designed to promote bone growth. TigerShark Threaded TL has a smooth bullet nose designed for fast insertion. TigerShark Threaded TL system must be used with supplemental fixation for the lumbar spine.

Features:

- BioBond™ porous structure technology designed to promote bone growth.
- Large open graft window to pack Stratofuse® biologics
- Simple threaded inserter
- Various sizes to accommodate different patient anatomies

Lordosis	Length	Width	Height
6°	25, 28, 32, 36mm	11mm	7-14mm (1mm increments)



400 Erin Drive, Knoxville, TN 37919 | O: 865.246.3333 | F: 865.246.3334 | choicespine.com